

For Research Use Only

# Fas/CD95 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00430-3

## Capture Antibody Information

Catalog Number:  
83400-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
240340E9  
Reactivity:  
human

Conjugate:  
Unconjugated  
Full name:  
Fas (TNF receptor superfamily,  
member 6)  
Gene ID:  
355

## Detection Antibody Information

Catalog Number:  
83400-1-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
240340B6  
Reactivity:  
human  
GenBank:  
NM\_000043.6

Conjugate:  
Unconjugated  
Full name:  
Fas (TNF receptor superfamily,  
member 6)  
Gene ID:  
355

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00430-3 targets Fas/CD95 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Fas/CD95 Recombinant antibody, PBS Only (Capture) 83400-2-PBS (240340E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: Fas/CD95 Recombinant antibody, PBS Only (Detector) 83400-1-PBS (240340B6). 100 µg. Concentration 1 mg/mL.

Alternative Fas/CD95 matched antibody pairs: MP00430-1, MP00430-2, MP50205-1, MP50205-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00430-3, Fas/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83400-2-PBS. Detection antibody: 83400-1-PBS. Standard: Eg0877. Range: 1.25-40 ng/mL