

For Research Use Only



# TSSK4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00436-2

## Capture Antibody Information

Catalog Number:  
83401-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240373D8

Reactivity:  
Human

Immunogen Catalog Number:  
Ag16205

Conjugate:  
Unconjugated

Full name:  
testis-specific serine kinase 4

Gene ID:  
283629

## Detection Antibody Information

Catalog Number:  
83401-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240373F10

Reactivity:  
Human

GenBank:  
BC111088

Immunogen Catalog Number:  
Ag16205

Conjugate:  
Unconjugated

Full name:  
testis-specific serine kinase 4

Gene ID:  
283629

## Applications

Tested Applications:  
Cytometric bead array

Range:  
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00436-2 targets TSSK4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSSK4 Recombinant antibody, PBS Only (Capture/Detector) 83401-1-PBS (240373D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: TSSK4 Recombinant antibody, PBS Only (Detector) 83401-4-PBS (240373F10). 100 µg. Concentration 1 mg/mL.

Alternative TSSK4 matched antibody pairs: MP00436-1, MP00436-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00436-2, TSSK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83401-1-PBS. Detection antibody: 83401-4-PBS. Standard: Ag16205. Range: 6.25-100 ng/mL