

TMEM119 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00444-3**

Capture Antibody Information

Catalog Number:
83448-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240433B9
Reactivity:
human

Conjugate:
Unconjugated
Full name:
transmembrane protein 119
Gene ID:
338773

Detection Antibody Information

Catalog Number:
83448-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240433C12
Reactivity:
human
GenBank:
NM_181724.2

Conjugate:
Unconjugated
Full name:
transmembrane protein 119
Gene ID:
338773

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00444-3 targets TMEM119 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMEM119 Recombinant antibody, PBS Only (Capture) 83448-1-PBS (240433B9). 100 µg. Concentration 1 mg/mL.

Detection antibody: TMEM119 Recombinant antibody, PBS Only (Capture/Detector) 83448-4-PBS (240433C12). 100 µg. Concentration 1 mg/mL.

Alternative TMEM119 matched antibody pairs: MP00444-1, MP00444-2, MP00444-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00444-3, TMEM119 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83448-1-PBS. Detection antibody: 83448-4-PBS. Standard: Eg0796. Range: 1.25-40 ng/mL