

# RBMS1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00575-1**

## Capture Antibody Information

**Catalog Number:**  
83623-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240677G10  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag1493

**Conjugate:**  
Unconjugated  
**Full name:**  
RNA binding motif, single stranded interacting protein 1  
**Gene ID:**  
5937

## Detection Antibody Information

**Catalog Number:**  
83623-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240677E5  
**Reactivity:**  
human  
**GenBank:**  
BC018951  
**Immunogen Catalog Number:**  
Ag1493

**Conjugate:**  
Unconjugated  
**Full name:**  
RNA binding motif, single stranded interacting protein 1  
**Gene ID:**  
5937

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.78-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00575-1 targets RBMS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBMS1 Recombinant antibody, PBS Only (Capture) 83623-4-PBS (240677G10). 100 µg. Concentration 1 mg/mL.

Detection antibody: RBMS1 Recombinant antibody, PBS Only (Detector) 83623-1-PBS (240677E5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00575-1, RBMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83623-4-PBS. Detection antibody: 83623-1-PBS. Standard: Ag1493. Range: 0.78-100 ng/mL