

# Mouse RAGE Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00693-2**

## Capture Antibody Information

**Catalog Number:**  
83742-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240774H4  
**Reactivity:**  
mouse  
**Immunogen Catalog Number:**  
Eg0873

**Conjugate:**  
Unconjugated  
**Full name:**  
advanced glycosylation end product-specific receptor  
**Gene ID:**  
11596

## Detection Antibody Information

**Catalog Number:**  
83742-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240774A12  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_007425.3  
**Immunogen Catalog Number:**  
Eg0873

**Conjugate:**  
Unconjugated  
**Full name:**  
advanced glycosylation end product-specific receptor  
**Gene ID:**  
11596

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.78-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00693-2 targets RAGE in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse RAGE Recombinant antibody, PBS Only (Detector) 83742-4-PBS (240774H4). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse RAGE Recombinant antibody, PBS Only 83742-1-PBS (240774A12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00693-2, Mouse RAGE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83742-4-PBS. Detection antibody: 83742-1-PBS. Standard: Eg0873. Range: 0.78-100 ng/mL.