

# Mesothelin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00697-1**

## Capture Antibody Information

Catalog Number:

83744-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240795C5

Reactivity:

human

Immunogen Catalog Number:

Eg31658

Conjugate:

Unconjugated

Full name:

mesothelin

Gene ID:

10232

## Detection Antibody Information

Catalog Number:

83744-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240795D7

Reactivity:

human

GenBank:

NM\_005823

Immunogen Catalog Number:

Eg31658

Conjugate:

Unconjugated

Full name:

mesothelin

Gene ID:

10232

## Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00697-1 targets Mesothelin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mesothelin Recombinant antibody, PBS Only (Capture) 83744-1-PBS (240795C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mesothelin Recombinant antibody, PBS Only (Detector) 83744-4-PBS (240795D7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00697-1, MSLN/Mesothelin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83744-1-PBS. Detection antibody: 83744-4-PBS. Standard: Eg31658. Range: 0.625-80 ng/mL.