

VEGF-C Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00710-2**

Capture Antibody Information

Catalog Number:
83727-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240772C1
Reactivity:
human

Conjugate:
Unconjugated
Full name:
vascular endothelial growth factor C
Gene ID:
7424

Detection Antibody Information

Catalog Number:
83727-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240772A5
Reactivity:
human
GenBank:
BC035212

Conjugate:
Unconjugated
Full name:
vascular endothelial growth factor C
Gene ID:
7424

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00710-2 targets VEGF-C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VEGF-C Recombinant antibody, PBS Only (Capture/Detector) 83727-3-PBS (240772C1). 100 µg. Concentration 1 mg/mL.

Detection antibody: VEGF-C Recombinant antibody, PBS Only (Detector) 83727-1-PBS (240772A5). 100 µg. Concentration 1 mg/mL.

Alternative VEGF-C matched antibody pairs: MP00710-1, MP00710-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00710-2, VEGF-C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83727-3-PBS. Detection antibody: 83727-1-PBS. Standard: Eg32025. Range: 0.313-20 ng/mL