

# HIF2 $\alpha$ /EPAS1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00775-1**

## Capture Antibody Information

**Catalog Number:**

83790-4-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240736H7

**Reactivity:**

human

**Immunogen Catalog Number:**

Ag15199

**Conjugate:**

Unconjugated

**Full name:**

endothelial PAS domain protein 1

**Gene ID:**

2034

## Detection Antibody Information

**Catalog Number:**

83790-1-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240736F2

**Reactivity:**

human, rat

**GenBank:**

BC051338

**Immunogen Catalog Number:**

Ag15199

**Conjugate:**

Unconjugated

**Full name:**

endothelial PAS domain protein 1

**Gene ID:**

2034

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.78-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00775-1 targets HIF2 $\alpha$ /EPAS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HIF2 $\alpha$ /EPAS1 Recombinant antibody, PBS Only (Capture) 83790-4-PBS (240736H7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: HIF2 $\alpha$ /EPAS1 Recombinant antibody, PBS Only (Detector) 83790-1-PBS (240736F2). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00775-1, HIF2α/EPAS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83790-4-PBS. Detection antibody: 83790-1-PBS. Standard: Ag15199. Range: 0.78-100 ng/mL.