

# FRMD3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00787-3**

## Capture Antibody Information

**Catalog Number:**

83850-1-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240967D5

**Reactivity:**

human

**Immunogen Catalog Number:**

Ag13235

**Conjugate:**

Unconjugated

**Full name:**

FERM domain containing 3

**Gene ID:**

257019

## Detection Antibody Information

**Catalog Number:**

83850-2-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240967D9

**Reactivity:**

human

**GenBank:**

BC023560

**Immunogen Catalog Number:**

Ag13235

**Conjugate:**

Unconjugated

**Full name:**

FERM domain containing 3

**Gene ID:**

257019

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.78-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00787-3 targets FRMD3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FRMD3 Recombinant antibody, PBS Only (Capture/Detector) 83850-1-PBS (240967D5). 100 µg. Concentration 1 mg/mL.

Detection antibody: FRMD3 Recombinant antibody, PBS Only (Detector) 83850-2-PBS (240967D9). 100 µg. Concentration 1 mg/mL.

Alternative FRMD3 matched antibody pairs: MP00787-1, MP00787-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00787-3, FRMD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83850-1-PBS. Detection antibody: 83850-2-PBS. Standard: Ag13235. Range: 0.78-100 ng/mL