

For Research Use Only

# Mouse Tissue Factor Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00859-2

## Capture Antibody Information

Catalog Number:  
83887-3-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
240875H1  
Reactivity:  
mouse

Conjugate:  
Unconjugated  
Full name:  
coagulation factor III  
Gene ID:  
14066

## Detection Antibody Information

Catalog Number:  
83887-4-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
240875H3  
Reactivity:  
human, mouse, rat  
GenBank:

Conjugate:  
Unconjugated  
Full name:  
coagulation factor III  
Gene ID:  
14066

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00859-2 targets Tissue Factor in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: Mouse Coagulation Factor III/Tissue Factor Recombinant antibody, PBS Only (Capture) 83887-3-PBS (240875H1). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Tissue Factor Recombinant antibody, PBS Only (Capture/Detector) 83887-4-PBS (240875H3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00859-2, MOUSE Tissue Factor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83887-3-PBS. Detection antibody: 83887-4-PBS. Standard: Eg1102. Range: 0.313-20 ng/mL