

HORMAD1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00934-3**

Capture Antibody Information

Catalog Number:
84003-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241117G1
Reactivity:
human
Immunogen Catalog Number:
Ag4914

Conjugate:
Unconjugated
Full name:
HORMA domain containing 1
Gene ID:
84072

Detection Antibody Information

Catalog Number:
84003-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241117A8
Reactivity:
human
GenBank:
BC047406
Immunogen Catalog Number:
Ag4914

Conjugate:
Unconjugated
Full name:
HORMA domain containing 1
Gene ID:
84072

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00934-3 targets HORMAD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HORMAD1 Recombinant antibody, PBS Only (Capture) 84003-4-PBS (241117G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: HORMAD1 Recombinant antibody, PBS Only (Detector) 84003-1-PBS (241117A8). 100 µg. Concentration 1 mg/mL.

Alternative HORMAD1 matched antibody pairs: MP00934-1, MP00934-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00934-3, HORMAD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84003-4-PBS. Detection antibody: 84003-1-PBS. Standard: Ag4914. Range: 0.313-40 ng/mL