

# CD36 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00992-1**

## Capture Antibody Information

**Catalog Number:**  
84078-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241262F8  
**Reactivity:**  
human

**Conjugate:**  
Unconjugated  
**Full name:**  
CD36 molecule (thrombospondin receptor)  
**Gene ID:**  
948

## Detection Antibody Information

**Catalog Number:**  
84078-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241262C1  
**Reactivity:**  
human  
**GenBank:**  
NM\_001001548.3

**Conjugate:**  
Unconjugated  
**Full name:**  
CD36 molecule (thrombospondin receptor)  
**Gene ID:**  
948

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.563-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00992-1 targets CD36 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD36 Recombinant antibody, PBS Only (Capture/Detector) 84078-3-PBS (241262F8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD36 Recombinant antibody, PBS Only (Detector) 84078-1-PBS (241262C1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00992-1, CD36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84078-3-PBS. Detection antibody: 84078-1-PBS. Standard: Eg1070. Range: 1.563-100 ng/mL.