

EPCAM/CD326 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01009-2**

Capture Antibody Information

Catalog Number:
84073-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241243F9
Reactivity:
human

Conjugate:
Unconjugated
Full name:
epithelial cell adhesion molecule
Gene ID:
4072

Detection Antibody Information

Catalog Number:
84073-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241243E4
Reactivity:
human
GenBank:
NM_002354.3

Conjugate:
Unconjugated
Full name:
epithelial cell adhesion molecule
Gene ID:
4072

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01009-2 targets EPCAM/CD326 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EPCAM/CD326 Recombinant antibody, PBS Only (Capture) 84073-5-PBS (241243F9). 100 µg. Concentration 1 mg/mL.

Detection antibody: EPCAM/CD326 Recombinant antibody, PBS Only (Detector) 84073-4-PBS (241243E4). 100 µg. Concentration 1 mg/mL.

Alternative EPCAM/CD326 matched antibody pairs: MP01009-1, MP01009-3, MP01009-4, MP01344-1, MP01344-2, MP01344-3, MP50309-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01009-2, EPCAM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84073-5-PBS. Detection antibody: 84073-4-PBS. Standard: Eg1370. Range: 0.156-20 ng/mL