

For Research Use Only



ACMSD Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01022-1

Capture Antibody Information

Catalog Number:
84100-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241077F1
Reactivity:
human
Immunogen Catalog Number:
Ag25212

Conjugate:
Unconjugated
Full name:
aminocarboxymuconate
semialdehyde decarboxylase
Gene ID:
130013

Detection Antibody Information

Catalog Number:
84100-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241077D7
Reactivity:
human
GenBank:
BC016018
Immunogen Catalog Number:
Ag25212

Conjugate:
Unconjugated
Full name:
aminocarboxymuconate
semialdehyde decarboxylase
Gene ID:
130013

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-80 ng/mL (Cytometric Bead Array)
Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01022-1 targets ACMSD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ACMSD Recombinant antibody, PBS Only (Capture) 84100-2-PBS (241077F1). 100 µg. Concentration 1 mg/mL.

Detection antibody: ACMSD Recombinant antibody, PBS Only (Detector) 84100-1-PBS (241077D7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01022-1, ACMSD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84100-2-PBS. Detection antibody: 84100-1-PBS. Standard: Ag25212. Range: 2.5-80 ng/mL.