

For Research Use Only



ZNF350 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01034-1

Capture Antibody Information

Catalog Number: 84121-4-PBS	Clone ID: 241225G1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: zinc finger protein 350
Isotype: IgG	Immunogen Catalog Number: Ag23982	Gene ID: 59348
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 84121-2-PBS	Clone ID: 241225B9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: zinc finger protein 350
Isotype: IgG	GenBank: BC009921	Gene ID: 59348
Purification Method: Protein A purification	Immunogen Catalog Number: Ag23982	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
-----------------------------------------------	---------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP01034-1 targets ZNF350 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZNF350 Recombinant antibody, PBS Only (Capture) 84121-4-PBS (241225G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: ZNF350 Recombinant antibody, PBS Only (Detector) 84121-2-PBS (241225B9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01034-1, ZNF350 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84121-4-PBS. Detection antibody: 84121-2-PBS. Standard: Ag23982. Range: 0.781-100 ng/mL.