

For Research Use Only



TMIGD1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01037-2

Capture Antibody Information

Catalog Number:
84130-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241399E8
Reactivity:
human
Immunogen Catalog Number:
Ag25691

Conjugate:
Unconjugated
Full name:
transmembrane and immunoglobulin domain containing 1
Gene ID:
388364

Detection Antibody Information

Catalog Number:
84130-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241399G9
Reactivity:
human
GenBank:
BC137201
Immunogen Catalog Number:
Ag25691

Conjugate:
Unconjugated
Full name:
transmembrane and immunoglobulin domain containing 1
Gene ID:
388364

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01037-2 targets TMIGD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMIGD1 Recombinant antibody, PBS Only (Capture) 84130-1-PBS (241399E8). 100 µg. Concentration 1 mg/mL.

Detection antibody: TMIGD1 Recombinant antibody, PBS Only (Capture/Detector) 84130-4-PBS (241399G9). 100 µg. Concentration 1 mg/mL.

Alternative TMIGD1 matched antibody pairs: MP01037-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01037-2, TMIGD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84130-1-PBS. Detection antibody: 84130-4-PBS. Standard: Ag25691. Range: 1.563-200 ng/mL