

# IL-23R Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01048-2**

## Capture Antibody Information

**Catalog Number:**  
84112-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240905C9  
**Reactivity:**  
human

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 23 receptor  
**Gene ID:**  
149233

## Detection Antibody Information

**Catalog Number:**  
84112-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240905D5  
**Reactivity:**  
human  
**GenBank:**  
NM\_144701.3

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 23 receptor  
**Gene ID:**  
149233

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01048-2 targets IL-23R in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL23R Recombinant antibody, PBS Only (Capture) 84112-1-PBS (240905C9). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL23R Recombinant antibody, PBS Only (Detector) 84112-2-PBS (240905D5). 100 µg. Concentration 1 mg/mL.

Alternative IL-23R matched antibody pairs: MP01048-1, MP01048-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01048-2, IL23R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84112-1-PBS. Detection antibody: 84112-2-PBS. Standard: Eg0853. Range: 0.625-80 ng/mL