

Mouse CD80 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01051-1**

Capture Antibody Information

Catalog Number:
84137-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241371C5
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
CD80 antigen
Gene ID:
12519

Detection Antibody Information

Catalog Number:
84137-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241371E12
Reactivity:
mouse
GenBank:
NM_009855.2

Conjugate:
Unconjugated
Full name:
CD80 antigen
Gene ID:
12519

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01051-1 targets CD80 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD80 Recombinant antibody, PBS Only (Capture) 84137-2-PBS (241371C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD80 Recombinant antibody, PBS Only (Capture/Detector) 84137-4-PBS (241371E12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01051-1, MOUSE CD80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84137-2-PBS. Detection antibody: 84137-4-PBS. Standard: Eg0840. Range: 0.313-40 ng/mL.