

RENBP Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01073-1**

Capture Antibody Information

Catalog Number: 84152-2-PBS	Clone ID: 241142E6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: renin binding protein
Isotype: IgG	Immunogen Catalog Number: Ag35056	Gene ID: 5973
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 84152-1-PBS	Clone ID: 241142A5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: renin binding protein
Isotype: IgG	GenBank: BC015558	Gene ID: 5973
Purification Method: Protein A purification	Immunogen Catalog Number: Ag35056	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP01073-1 targets RENBP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RENBP Recombinant antibody, PBS Only (Capture) 84152-2-PBS (241142E6). 100 µg. Concentration 1 mg/mL.

Detection antibody: RENBP Recombinant antibody, PBS Only (Detector) 84152-1-PBS (241142A5). 100 µg. Concentration 1 mg/mL.

Alternative RENBP matched antibody pairs: MP01073-2, MP01073-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01073-1, RENBP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84152-2-PBS. Detection antibody: 84152-1-PBS. Standard: Ag35056. Range: 0.781-100 ng/mL.