

For Research Use Only

# TIM3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01083-3

## Capture Antibody Information

Catalog Number:  
84183-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
241421A10  
Reactivity:  
human

Conjugate:  
Unconjugated  
Full name:  
hepatitis A virus cellular receptor 2  
Gene ID:  
84868

## Detection Antibody Information

Catalog Number:  
84183-1-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
241421A3  
Reactivity:  
human  
GenBank:  
BC020843

Conjugate:  
Unconjugated  
Full name:  
hepatitis A virus cellular receptor 2  
Gene ID:  
84868

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01083-3 targets TIM3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HAVCR2 Recombinant antibody, PBS Only (Capture) 84183-2-PBS (241421A10). 100 µg. Concentration 1 mg/mL.

Detection antibody: HAVCR2 Recombinant antibody, PBS Only (Detector) 84183-1-PBS (241421A3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01083-3, TIM3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84183-2-PBS. Detection antibody: 84183-1-PBS. Standard: RP02340. Range: 0.781-100 ng/mL.