

# EPS15 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01095-3**

## Capture Antibody Information

**Catalog Number:**  
84179-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241485G9  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag35282

**Conjugate:**  
Unconjugated  
**Full name:**  
epidermal growth factor receptor pathway substrate 15  
**Gene ID:**  
2060

## Detection Antibody Information

**Catalog Number:**  
84179-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241485G5  
**Reactivity:**  
human  
**GenBank:**  
NM\_001981.3  
**Immunogen Catalog Number:**  
Ag35282

**Conjugate:**  
Unconjugated  
**Full name:**  
epidermal growth factor receptor pathway substrate 15  
**Gene ID:**  
2060

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01095-3 targets EPS15 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EPS15 Recombinant antibody, PBS Only (Capture/Detector) 84179-4-PBS (241485G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: EPS15 Recombinant antibody, PBS Only (Detector) 84179-3-PBS (241485G5). 100 µg. Concentration 1 mg/mL.

Alternative EPS15 matched antibody pairs: MP01095-1, MP01095-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01095-3, EPS15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84179-4-PBS. Detection antibody: 84179-3-PBS. Standard: Ag35282. Range: 0.625-80 ng/mL