

PYROXD1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01201-2**

Capture Antibody Information

Catalog Number:
84289-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241639C6
Reactivity:
human
Immunogen Catalog Number:
Ag34768

Conjugate:
Unconjugated
Full name:
pyridine nucleotide-disulphide oxidoreductase domain 1
Gene ID:
79912

Detection Antibody Information

Catalog Number:
84289-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241639B6
Reactivity:
human
GenBank:
BC021662
Immunogen Catalog Number:
Ag34768

Conjugate:
Unconjugated
Full name:
pyridine nucleotide-disulphide oxidoreductase domain 1
Gene ID:
79912

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01201-2 targets PYROXD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PYROXD1 Recombinant antibody, PBS Only (Capture) 84289-3-PBS (241639C6). 100 µg. Concentration 1 mg/mL.

Detection antibody: PYROXD1 Recombinant antibody, PBS Only (Detector) 84289-1-PBS (241639B6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01201-2, PYROXD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84289-3-PBS. Detection antibody: 84289-1-PBS. Standard: Ag34768. Range: 0.313-40 ng/mL.