

CSH1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01252-3**

Capture Antibody Information

Catalog Number:
84361-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241625F7
Reactivity:
human
Immunogen Catalog Number:
Ag0201

Conjugate:
Unconjugated
Full name:
chorionic somatomammotropin hormone 1 (placental lactogen)
Gene ID:
1442

Detection Antibody Information

Catalog Number:
84361-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241625C9
Reactivity:
human
GenBank:
BC002717
Immunogen Catalog Number:
Ag0201

Conjugate:
Unconjugated
Full name:
chorionic somatomammotropin hormone 1 (placental lactogen)
Gene ID:
1442

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01252-3 targets CSH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Placental Lactogen Recombinant antibody, PBS Only (Capture) 84361-3-PBS (241625F7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Placental Lactogen Recombinant antibody, PBS Only (Detector) 84361-1-PBS (241625C9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01252-3, CSH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84361-3-PBS. Detection antibody: 84361-1-PBS. Standard: Ag0201. Range: 0.781-100 ng/mL.