

# Mouse Nectin-2/CD112

## Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01292-3**

### Capture Antibody Information

**Catalog Number:**  
84439-5-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241798F12  
**Reactivity:**  
mouse

**Conjugate:**  
Unconjugated  
**Full name:**  
poliovirus receptor-related 2  
**Gene ID:**  
19294

### Detection Antibody Information

**Catalog Number:**  
84439-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241798E6  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_008990.3

**Conjugate:**  
Unconjugated  
**Full name:**  
poliovirus receptor-related 2  
**Gene ID:**  
19294

### Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

### Product Information

MP01292-3 targets Nectin-2/CD112 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Nectin-2/CD112 Recombinant antibody, PBS Only (Capture) 84439-5-PBS (241798F12). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Nectin-2/CD112 Recombinant antibody, PBS Only (Detector) 84439-4-PBS (241798E6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

### Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01292-3, MOUSE Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84439-5-PBS. Detection antibody: 84439-4-PBS. Standard: Eg2019. Range: 0.313-20 ng/mL