

Rat IL-4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01293-3**

Capture Antibody Information

Catalog Number:
84422-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241249H7
Reactivity:
rat

Conjugate:
Unconjugated
Full name:
interleukin 4
Gene ID:
287287

Detection Antibody Information

Catalog Number:
84422-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241249D5
Reactivity:
rat
GenBank:
NM_201270.1

Conjugate:
Unconjugated
Full name:
interleukin 4
Gene ID:
287287

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01293-3 targets IL-4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat IL-4 Recombinant antibody, PBS Only (Capture/Detector) 84422-5-PBS (241249H7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Rat IL-4 Recombinant antibody, PBS Only (Detector) 84422-2-PBS (241249D5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01293-3, RAT IL-4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-5-PBS. Detection antibody: 84422-2-PBS. Standard: Eg1113. Range: 1.25-80 ng/mL