

For Research Use Only



POLR2K Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01314-1

Capture Antibody Information

Catalog Number:
82995-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230147G10

Reactivity:
human

Immunogen Catalog Number:
Ag33765

Conjugate:
Unconjugated

Full name:
polymerase (RNA) II (DNA directed)
polypeptide K, 7.0kDa

Gene ID:
5440

Detection Antibody Information

Catalog Number:
82995-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230147E8

Reactivity:
human

GenBank:
BC018157

Immunogen Catalog Number:
Ag33765

Conjugate:
Unconjugated

Full name:
polymerase (RNA) II (DNA directed)
polypeptide K, 7.0kDa

Gene ID:
5440

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01314-1 targets POLR2K in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: POLR2K Recombinant antibody, PBS Only (Capture) 82995-3-PBS (230147G10). 100 µg. Concentration 1 mg/mL.

Detection antibody: POLR2K Recombinant antibody, PBS Only (Detector) 82995-1-PBS (230147E8). 100 µg. Concentration 1 mg/mL.

Alternative POLR2K matched antibody pairs: MP01314-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01314-1, POLR2K Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82995-3-PBS. Detection antibody: 82995-1-PBS. Standard: Ag33765. Range: 0.313-20 ng/mL