

CLEC4A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01361-1**

Capture Antibody Information

Catalog Number:
84515-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241854G1
Reactivity:
human

Conjugate:
Unconjugated
Full name:
C-type lectin domain family 4, member A
Gene ID:
50856

Detection Antibody Information

Catalog Number:
84515-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241854B10
Reactivity:
human
GenBank:
NM_016184.4

Conjugate:
Unconjugated
Full name:
C-type lectin domain family 4, member A
Gene ID:
50856

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01361-1 targets CLEC4A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DCIR/CLEC4A Recombinant antibody, PBS Only (Capture) 84515-4-PBS (241854G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: DCIR/CLEC4A Recombinant antibody, PBS Only (Detector) 84515-1-PBS (241854B10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01361-1, CLEC4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84515-4-PBS. Detection antibody: 84515-1-PBS. Standard: Eg1350. Range: 0.469-60 ng/mL.