

# IL-2RG Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01367-1**

## Capture Antibody Information

**Catalog Number:**  
84501-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242020H6  
**Reactivity:**  
human

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 2 receptor, gamma  
(severe combined immunodeficiency)  
**Gene ID:**  
3561

## Detection Antibody Information

**Catalog Number:**  
84501-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242020G10  
**Reactivity:**  
human  
**GenBank:**  
NM\_000206.3

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 2 receptor, gamma  
(severe combined immunodeficiency)  
**Gene ID:**  
3561

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01367-1 targets IL-2RG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-2 R gamma/CD132 Recombinant antibody, PBS Only (Capture) 84501-2-PBS (242020H6). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-2 R gamma/CD132 Recombinant antibody, PBS Only (Detector) 84501-4-PBS (242020G10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01367-1, IL-2RG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84501-2-PBS. Detection antibody: 84501-4-PBS. Standard: Eg1705. Range: 0.313-40 ng/mL