

RAD51C Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01399-2**

Capture Antibody Information

Catalog Number:

84537-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241971D4

Reactivity:

human

Immunogen Catalog Number:

Ag35070

Conjugate:

Unconjugated

Full name:

RAD51 homolog C (*S. cerevisiae*)

Gene ID:

5889

Detection Antibody Information

Catalog Number:

84537-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241971F6

Reactivity:

human

GenBank:

BC107753

Immunogen Catalog Number:

Ag35070

Conjugate:

Unconjugated

Full name:

RAD51 homolog C (*S. cerevisiae*)

Gene ID:

5889

Applications

Tested Applications:

Cytometric bead array

Range:

0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01399-2 targets RAD51C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAD51C Recombinant antibody, PBS Only (Capture/Detector) 84537-1-PBS (241971D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: RAD51C Recombinant antibody, PBS Only (Detector) 84537-3-PBS (241971F6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01399-2, RAD51C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84537-1-PBS. Detection antibody: 84537-3-PBS. Standard: Ag35070. Range: 0.469-60 ng/mL.