

EMR2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01442-3**

Capture Antibody Information

Catalog Number:
84608-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242042B11
Reactivity:
human

Conjugate:
Unconjugated
Full name:
egf-like module containing, mucin-like, hormone receptor-like 2
Gene ID:
30817

Detection Antibody Information

Catalog Number:
84608-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242042G9
Reactivity:
human
GenBank:
NM_013447.3

Conjugate:
Unconjugated
Full name:
egf-like module containing, mucin-like, hormone receptor-like 2
Gene ID:
30817

Applications

Tested Applications:
Sandwich ELISA

Range:
39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01442-3 targets EMR2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: EMR2 Recombinant antibody, PBS Only (Capture) 84608-2-PBS (242042B11). 100 µg. Concentration 1 mg/mL.

Detection antibody: EMR2 Recombinant antibody, PBS Only (Detector) 84608-4-PBS (242042G9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01442-3, Human EMR2 Recombinant Matched Antibody Pair - PBS only. 84608-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2479. 84608-4-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.



The mean EMR2 concentration was determined to be 93.8 ng/mL in HL-60 cell extract based on a 1.4 mg/mL extract load, 98.6 ng/mL in THP-1 cell extract based on a 4.4 mg/mL extract load, 27.8 ng/mL in K-562 cell extract based on a 2.4 mg/mL extract load and 28.5 ng/mL in Jurkat cell extract based on a 4.2 mg/mL extract load.