

NEK4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01486-2**

Capture Antibody Information

Catalog Number:
84654-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242081A8
Reactivity:
human
Immunogen Catalog Number:
Ag34929

Conjugate:
Unconjugated
Full name:
NIMA (never in mitosis gene a)-related kinase 4
Gene ID:
6787

Detection Antibody Information

Catalog Number:
84654-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242081F4
Reactivity:
human
GenBank:
BC063044
Immunogen Catalog Number:
Ag34929

Conjugate:
Unconjugated
Full name:
NIMA (never in mitosis gene a)-related kinase 4
Gene ID:
6787

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01486-2 targets NEK4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NEK4 Recombinant antibody, PBS Only (Capture) 84654-1-PBS (242081A8). 100 µg. Concentration 1 mg/mL.

Detection antibody: NEK4 Recombinant antibody, PBS Only (Detector) 84654-3-PBS (242081F4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01486-2, NEK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84654-1-PBS. Detection antibody: 84654-3-PBS. Standard: Ag34929. Range: 0.625-80 ng/mL