

CYP51A1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01526-2**

Capture Antibody Information

Catalog Number:
84705-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242166D6
Reactivity:
human
Immunogen Catalog Number:
Ag4235

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 51, subfamily A, polypeptide 1
Gene ID:
1595

Detection Antibody Information

Catalog Number:
84705-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242166G8
Reactivity:
human
GenBank:
BC032322
Immunogen Catalog Number:
Ag4235

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 51, subfamily A, polypeptide 1
Gene ID:
1595

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01526-2 targets CYP51A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CYP51A1 Recombinant antibody, PBS Only (Capture) 84705-2-PBS (242166D6). 100 µg. Concentration 1 mg/mL.

Detection antibody: CYP51A1 Recombinant antibody, PBS Only (Detector) 84705-3-PBS (242166G8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01526-2, CYP51A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84705-2-PBS. Detection antibody: 84705-3-PBS. Standard: Ag4235. Range: 0.781-100 ng/mL.