

# TMIGD2/CD28H Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01528-1**

## Capture Antibody Information

**Catalog Number:**  
84727-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241807D1  
**Reactivity:**  
human

**Conjugate:**  
Unconjugated  
**Full name:**  
transmembrane and immunoglobulin  
domain containing 2  
**Gene ID:**  
126259

## Detection Antibody Information

**Catalog Number:**  
84727-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241807B12  
**Reactivity:**  
human  
**GenBank:**  
BC015655

**Conjugate:**  
Unconjugated  
**Full name:**  
transmembrane and immunoglobulin  
domain containing 2  
**Gene ID:**  
126259

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead  
Array)

**Recommended Dilutions:**  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP01528-1 targets TMIGD2/CD28H in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMIGD2/CD28H Recombinant antibody, PBS Only (Capture) 84727-4-PBS (241807D1). 100 µg. Concentration 1 mg/mL.

Detection antibody: TMIGD2/CD28H Recombinant antibody, PBS Only (Detector) 84727-2-PBS (241807B12). 100 µg. Concentration 1 mg/mL.

Alternative TMIGD2/CD28H matched antibody pairs: MP01528-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01528-1, TMIGD2/CD28H Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84727-4-PBS. Detection antibody: 84727-2-PBS. Standard: Eg2264. Range: 0.781-100 ng/mL.