

cIAP2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01556-2**

Capture Antibody Information

Catalog Number:
84779-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242214G1
Reactivity:
human
Immunogen Catalog Number:
Ag20500

Conjugate:
Unconjugated
Full name:
baculoviral IAP repeat-containing 3
Gene ID:
330

Detection Antibody Information

Catalog Number:
84779-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242214E12
Reactivity:
human
GenBank:
BC027485
Immunogen Catalog Number:
Ag20500

Conjugate:
Unconjugated
Full name:
baculoviral IAP repeat-containing 3
Gene ID:
330

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01556-2 targets cIAP2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: cIAP2 Recombinant antibody, PBS Only (Capture) 84779-3-PBS (242214G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: cIAP2 Recombinant antibody, PBS Only (Capture/Detector) 84779-2-PBS (242214E12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01556-2, BIRC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84779-3-PBS. Detection antibody: 84779-2-PBS. Standard: Ag20500. Range: 0.313-20 ng/mL