

# SPTAN1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01560-1**

## Capture Antibody Information

**Catalog Number:**  
84761-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242104B7  
**Reactivity:**  
human, mouse, rat  
**Immunogen Catalog Number:**  
Ag36472

**Conjugate:**  
Unconjugated  
**Full name:**  
spectrin, alpha, non-erythrocytic 1 (alpha-fodrin)  
**Gene ID:**  
6709

## Detection Antibody Information

**Catalog Number:**  
84761-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242104B11  
**Reactivity:**  
human  
**GenBank:**  
NM\_003127.3  
**Immunogen Catalog Number:**  
Ag36472

**Conjugate:**  
Unconjugated  
**Full name:**  
spectrin, alpha, non-erythrocytic 1 (alpha-fodrin)  
**Gene ID:**  
6709

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.469-60 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01560-1 targets SPTAN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPTAN1 Recombinant antibody, PBS Only 84761-2-PBS (242104B7). 100 µg. Concentration 1 mg/mL.

Detection antibody: SPTAN1 Recombinant antibody, PBS Only 84761-1-PBS (242104B11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01560-1, SPTAN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84761-2-PBS. Detection antibody: 84761-1-PBS. Standard: Ag36472. Range: 0.469-60 ng/mL.