

SLC3A1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01622-2**

Capture Antibody Information

Catalog Number:
84856-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242370D3
Reactivity:
human, mouse, rat
Immunogen Catalog Number:
Ag9525

Conjugate:
Unconjugated
Full name:
solute carrier family 3 (cystine, dibasic and neutral amino acid transporters, activator of cystine, dibasic and neutral amino acid transport), member 1
Gene ID:
6519

Detection Antibody Information

Catalog Number:
84856-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242370H1
Reactivity:
human
GenBank:
BC022386
Immunogen Catalog Number:
Ag9525

Conjugate:
Unconjugated
Full name:
solute carrier family 3 (cystine, dibasic and neutral amino acid transporters, activator of cystine, dibasic and neutral amino acid transport), member 1
Gene ID:
6519

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01622-2 targets SLC3A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC3A1 Recombinant antibody, PBS Only (Capture) 84856-2-PBS (242370D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC3A1 Recombinant antibody, PBS Only (Detector) 84856-3-PBS (242370H1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01622-2, SLC3A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84856-2-PBS. Detection antibody: 84856-3-PBS. Standard: Ag9525. Range: 0.625-80 ng/mL