

Mouse OX40 Ligand/TNFSF4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01856-1**

Capture Antibody Information

Catalog Number:
85105-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242227F1
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand)
superfamily, member 4
Gene ID:
22164

Detection Antibody Information

Catalog Number:
85105-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242227G2
Reactivity:
mouse
GenBank:
NM_009452.2

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand)
superfamily, member 4
Gene ID:
22164

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-200 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01856-1 targets OX40 Ligand/TNFSF4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse OX40 Ligand/TNFSF4 Recombinant antibody, PBS Only (Capture) 85105-1-PBS (242227F1). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse OX40 Ligand/TNFSF4 Recombinant antibody, PBS Only (Capture/Detector) 85105-2-PBS (242227G2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



Cytometric bead array standard curve of MP01856-1, MOUSE OX40 Ligand/TNFSF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85105-1-PBS. Detection antibody: 85105-2-PBS. Standard: Eg1573. Range: 3.125-200 ng/mL