

PCBP3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01894-2**

Capture Antibody Information

Catalog Number:

85180-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242791F11

Reactivity:

human, mouse, rat

Immunogen Catalog Number:

Ag9202

Conjugate:

Unconjugated

Full name:

poly(rC) binding protein 3

Gene ID:

54039

Detection Antibody Information

Catalog Number:

85180-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242791E12

Reactivity:

human

GenBank:

BC012061

Immunogen Catalog Number:

Ag9202

Conjugate:

Unconjugated

Full name:

poly(rC) binding protein 3

Gene ID:

54039

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01894-2 targets PCBP3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCBP3 Recombinant antibody, PBS Only (Capture/Detector) 85180-3-PBS (242791F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: PCBP3 Recombinant antibody, PBS Only (Detector) 85180-1-PBS (242791E12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01894-2, PCBP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85180-3-PBS. Detection antibody: 85180-1-PBS. Standard: Ag9202. Range: 0.781-100 ng/mL.