

DFFA/DFF45 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01914-1**

Capture Antibody Information

Catalog Number:
85220-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242915G9
Reactivity:
human
Immunogen Catalog Number:
Ag0244

Conjugate:
Unconjugated
Full name:
DNA fragmentation factor, 45kDa, alpha polypeptide
Gene ID:
1676

Detection Antibody Information

Catalog Number:
85220-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242915C7
Reactivity:
human
GenBank:
BC000037
Immunogen Catalog Number:
Ag0244

Conjugate:
Unconjugated
Full name:
DNA fragmentation factor, 45kDa, alpha polypeptide
Gene ID:
1676

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01914-1 targets DFFA/DFF45 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: DFFA,DFF45 Recombinant antibody, PBS Only (Capture) 85220-3-PBS (242915G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: DFFA,DFF45 Recombinant antibody, PBS Only (Detector) 85220-1-PBS (242915C7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01914-1, DFFA/DFF45 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85220-3-PBS. Detection antibody: 85220-1-PBS. Standard: Ag0244. Range: 0.781-100 ng/mL.