

DLL4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01916-1**

Capture Antibody Information

Catalog Number: 85216-1-PBS	Clone ID: 242833E3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: delta-like 4 (Drosophila)
Isotype: IgG	Immunogen Catalog Number: Ag16140	Gene ID: 54567
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 85216-3-PBS	Clone ID: 242833G3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: delta-like 4 (Drosophila)
Isotype: IgG	GenBank: BC106950	Gene ID: 54567
Purification Method: Protein A purification	Immunogen Catalog Number: Ag16140	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP01916-1 targets DLL4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DLL4 Recombinant antibody, PBS Only (Capture) 85216-1-PBS (242833E3). 100 µg. Concentration 1 mg/mL.

Detection antibody: DLL4 Recombinant antibody, PBS Only (Capture/Detector) 85216-3-PBS (242833G3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01916-1, DLL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85216-1-PBS. Detection antibody: 85216-3-PBS. Standard: Ag16140. Range: 0.781-100 ng/mL.