

For Research Use Only



SRPK2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01968-1

Capture Antibody Information

Catalog Number: 85613-3-PBS	Clone ID: 242935C10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: SFRS protein kinase 2
Isotype: IgG	Immunogen Catalog Number: Ag17885	Gene ID: 6733
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 85613-1-PBS	Clone ID: 242935C7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: SFRS protein kinase 2
Isotype: IgG	GenBank: BC035214	Gene ID: 6733
Purification Method: Protein A purification	Immunogen Catalog Number: Ag17885	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP01968-1 targets SRPK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SRPK2 Recombinant antibody, PBS Only (Capture) 85613-3-PBS (242935C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SRPK2 Recombinant antibody, PBS Only (Detector) 85613-1-PBS (242935C7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01968-1, SRPK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85613-3-PBS. Detection antibody: 85613-1-PBS. Standard: Ag17885. Range: 0.781-100 ng/mL.