

# JAM2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01977-2**

## Capture Antibody Information

**Catalog Number:**

85601-4-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

242617B5

**Reactivity:**

human

**Conjugate:**

Unconjugated

**Full name:**

junctional adhesion molecule 2

**Gene ID:**

58494

## Detection Antibody Information

**Catalog Number:**

85601-3-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

242617H2

**Reactivity:**

human

**GenBank:**

NM\_021219.3

**Conjugate:**

Unconjugated

**Full name:**

junctional adhesion molecule 2

**Gene ID:**

58494

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01977-2 targets JAM2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: JAM2 Recombinant antibody, PBS Only (Capture) 85601-4-PBS (242617B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: JAM2 Recombinant antibody, PBS Only (Detector) 85601-3-PBS (242617H2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01977-2, JAM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85601-4-PBS. Detection antibody: 85601-3-PBS. Standard: Eg2914. Range: 0.313-40 ng/mL