

Mouse NGFR/TNFRSF16

Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02004-2**

Capture Antibody Information

Catalog Number:
85622-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242634F10
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
nerve growth factor receptor (TNFR
superfamily, member 16)
Gene ID:
18053

Detection Antibody Information

Catalog Number:
85622-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242634C9
Reactivity:
mouse
GenBank:
NM_033217.3

Conjugate:
Unconjugated
Full name:
nerve growth factor receptor (TNFR
superfamily, member 16)
Gene ID:
18053

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP02004-2 targets NGFR/TNFRSF16 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Ngfr Recombinant antibody, PBS Only (Capture/Detector) 85622-2-PBS (242634F10). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Ngfr Recombinant antibody, PBS Only (Detector) 85622-1-PBS (242634C9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



Cytometric bead array standard curve of MP02004-2, MOUSE NGFR/TNFRSF16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85622-2-PBS. Detection antibody: 85622-1-PBS. Standard: Eg2651. Range: 0.781-100 ng/mL.