

# SKP2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02007-3**

## Capture Antibody Information

**Catalog Number:**  
85670-7-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
243080C10  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag6826

**Conjugate:**  
Unconjugated  
**Full name:**  
S-phase kinase-associated protein 2 (p45)  
**Gene ID:**  
6502

## Detection Antibody Information

**Catalog Number:**  
85670-6-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
243080C11  
**Reactivity:**  
human  
**GenBank:**  
BC001441  
**Immunogen Catalog Number:**  
Ag6826

**Conjugate:**  
Unconjugated  
**Full name:**  
S-phase kinase-associated protein 2 (p45)  
**Gene ID:**  
6502

## Applications

**Tested Applications:**  
Sandwich ELISA

**Range:**  
0.625-40 ng/mL (Sandwich ELISA)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02007-3 targets SKP2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: SKP2 Recombinant antibody, PBS Only (Capture) 85670-7-PBS (243080C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SKP2 Recombinant antibody, PBS Only (Detector) 85670-6-PBS (243080C11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP02007-3, Human SKP2 Recombinant Matched Antibody Pair - PBS only. 85670-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag6826. 85670-6-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL



The mean SKP2 concentration was determined to be 21.92 ng/mL in U2OS cell extract based on a 1.40 mg/mL extract load and 29.08 ng/mL in HeLa cell extract based on a 1.80 mg/mL extract load.