

TCEA1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02021-1**

Capture Antibody Information

Catalog Number:
85664-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243155A9
Reactivity:
human
Immunogen Catalog Number:
Ag12199

Conjugate:
Unconjugated
Full name:
transcription elongation factor A (SII), 1
Gene ID:
6917

Detection Antibody Information

Catalog Number:
85664-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243155C9
Reactivity:
human
GenBank:
BC072460
Immunogen Catalog Number:
Ag12199

Conjugate:
Unconjugated
Full name:
transcription elongation factor A (SII), 1
Gene ID:
6917

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02021-1 targets TCEA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SII-TFIIS, TCEA1 Recombinant antibody, PBS Only (Capture) 85664-1-PBS (243155A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: SII-TFIIS, TCEA1 Recombinant antibody, PBS Only (Capture/Detector) 85664-3-PBS (243155C9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02021-1, TCEA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85664-1-PBS. Detection antibody: 85664-3-PBS. Standard: Ag12199. Range: 0.781-50 ng/mL