

ZEB1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02086-1**

Capture Antibody Information

Catalog Number:

85729-6-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

243139C5

Reactivity:

human

Immunogen Catalog Number:

Ag16204

Conjugate:

Unconjugated

Full name:

zinc finger E-box binding homeobox 1

Gene ID:

6935

Detection Antibody Information

Catalog Number:

85729-5-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

243139H4

Reactivity:

human

GenBank:

BC112392

Immunogen Catalog Number:

Ag16204

Conjugate:

Unconjugated

Full name:

zinc finger E-box binding homeobox 1

Gene ID:

6935

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA

Range:

0.625-80 ng/mL (Cytometric Bead Array)
39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02086-1 targets ZEB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: ZEB1 Recombinant antibody, PBS Only (Capture) 85729-6-PBS (243139C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: ZEB1 Recombinant antibody, PBS Only (Detector) 85729-5-PBS (243139H4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP02086-1, Human ZEB1 Recombinant Matched Antibody Pair - PBS only. 85729-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag16204. 85729-5-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



The mean ZEB1 concentration was determined to be 18.65 ng/mL in Jurkat cell extract based on a 1.9 mg/mL extract load and 2.41 ng/mL in MDA-MB-453 cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02086-1, ZEB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85729-6-PBS. Detection antibody: 85729-5-PBS. Standard: Ag16204. Range: 0.625-80 ng/mL