

Mouse Tnfsf18 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02134-2**

Capture Antibody Information

Catalog Number:
85842-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242937G3
Reactivity:
mouse
Immunogen Catalog Number:
Eg2850

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand)
superfamily, member 18
Gene ID:
240873

Detection Antibody Information

Catalog Number:
85842-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242937C4
Reactivity:
mouse
GenBank:
NM_183391.3
Immunogen Catalog Number:
Eg2850

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand)
superfamily, member 18
Gene ID:
240873

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02134-2 targets Tnfsf18 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Tnfsf18 Recombinant antibody, PBS Only (Capture) 85842-3-PBS (242937G3). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Tnfsf18 Recombinant antibody, PBS Only 85842-1-PBS (242937C4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02134-2, MOUSE Tnfsf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85842-3-PBS. Detection antibody: 85842-1-PBS. Standard: Eg2850. Range: 0.781-25 ng/mL.