

PLEKHH2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02268-1**

Capture Antibody Information

Catalog Number:
86166-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250291C2
Reactivity:
human
Immunogen Catalog Number:
Ag5409

Conjugate:
Unconjugated
Full name:
pleckstrin homology domain containing, family H (with MyTH4 domain) member 2
Gene ID:
130271

Detection Antibody Information

Catalog Number:
86166-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250291E10
Reactivity:
human
GenBank:
BC063310
Immunogen Catalog Number:
Ag5409

Conjugate:
Unconjugated
Full name:
pleckstrin homology domain containing, family H (with MyTH4 domain) member 2
Gene ID:
130271

Applications

Tested Applications:
Sandwich ELISA

Range:
0.625-20 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02268-1 targets PLEKHH2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: PLEKHH2 Recombinant antibody, PBS Only (Capture) 86166-2-PBS (250291C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: PLEKHH2 Recombinant antibody, PBS Only (Detector) 86166-1-PBS (250291E10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02268-1, Human PLEKHH2 Recombinant Matched Antibody Pair - PBS only. 86166-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5409. 86166-1-PBS was HRP conjugated as the detection antibody. Range: 0.625-20 ng/mL



The mean PLEKHH2 concentration was determined to be 5.99 ng/mL in Jurkat cell extract based on a 1.90 mg/mL extract load.