

For Research Use Only



# ANGPTL3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP02294-1

## Capture Antibody Information

Catalog Number:  
86193-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
250808A9

Reactivity:  
human

Immunogen Catalog Number:  
Ag2573

Conjugate:  
Unconjugated

Full name:  
angiopoietin-like 3

Gene ID:  
27329

## Detection Antibody Information

Catalog Number:  
86193-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
250808F7

Reactivity:  
human

GenBank:  
BC058287

Immunogen Catalog Number:  
Ag2573

Conjugate:  
Unconjugated

Full name:  
angiopoietin-like 3

Gene ID:  
27329

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02294-1 targets ANGPTL3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ANGPTL3 Recombinant antibody, PBS Only (Capture) 86193-2-PBS (250808A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: ANGPTL3 Recombinant antibody, PBS Only (Detector) 86193-1-PBS (250808F7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean ANGPTL3 concentration was determined to be 14.8 ng/mL in human serum.



Cytometric bead array standard curve of MP02294-1, ANGPTL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86193-2-PBS. Detection antibody: 86193-1-PBS. Standard: Ag2573. Range: 1.563-200 ng/mL.