

# FKBP1A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02313-1**

## Capture Antibody Information

**Catalog Number:**

86249-3-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

250923D6

**Reactivity:**

human

**Immunogen Catalog Number:**

Ag0406

**Conjugate:**

Unconjugated

**Full name:**

FK506 binding protein 1A, 12kDa

**Gene ID:**

2280

## Detection Antibody Information

**Catalog Number:**

86249-2-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

250923G1

**Reactivity:**

human

**GenBank:**

BC001925

**Immunogen Catalog Number:**

Ag0406

**Conjugate:**

Unconjugated

**Full name:**

FK506 binding protein 1A, 12kDa

**Gene ID:**

2280

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02313-1 targets FKBP1A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FKBP1A Recombinant antibody, PBS Only (Capture) 86249-3-PBS (250923D6). 100 µg. Concentration 1 mg/mL.

Detection antibody: FKBP1A Recombinant antibody, PBS Only (Detector) 86249-2-PBS (250923G1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean FKBP1A concentration was determined to be 149.2 ng/mL in Jurkat cell extract based on a 2.2 mg/mL extract load.



Cytometric bead array standard curve of MP02313-1, FKBP1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86249-3-PBS. Detection antibody: 86249-2-PBS. Standard: Ag0406. Range: 0.781-100 ng/mL.